



Name: \_\_\_\_\_

Class: \_\_\_\_\_

Date: \_\_\_\_\_



Earth Sciences

High School

9 hours

## Objective

To explore wind and hydrogen fuel cell power sources and try to improve the power output of both.

## Materials

- Horizon's Wind to Hydrogen Science Kit
- Electric Fan
- Distilled Water
- Stopwatch
- Horizon Renewable Energy Monitor (optional)

## Background

Humans have been using the wind to power machines since ancient times. Even today, you can see windmills powering water pumps on farms and the huge wind turbines that are becoming an important part of our electricity production. Read more about how humans have used wind energy in [History of Wind Power](#).

Wind turbines can be classified as Horizontal axis or vertical axis, depending on how they turn. Horizontal axis wind turbines (HAWTs), like the ones pictured in Figure 1, consist of a rotor and a horizontal main rotor at the top of a tower. They also typically have a gearbox, which turns a slower rotation into a faster one, and a generator which produces electricity. You can find out more about how wind turbines work by reading [Parts of a Wind Turbine](#).



HAWT pumping water on a

## farm

The blades of a turbine rely on the same principles as a plane's wing, with lift and drag affecting how they move. The forces involved in spinning turbine blades are discussed in detail in [Aerodynamics of Wind Turbines](#).



HAWT for electricity generation

Our wind turbine will also be used to produce hydrogen from water with an electrolyzer, which is also a fuel cell. You can find out how our electrolyzer works by reading [Introduction to Electrolyzers and Fuel Cells](#).

Wind turbines today are one of the ways we are harnessing clean energy to generate electricity. They create no pollution, but they are limited by where they can be used and when the wind is blowing. With these restrictions, a wind turbine must be built to generate the most electricity possible.

How can we make the wind turbine more effective at generating power? How does the power output of a wind turbine compare to a hydrogen fuel cell? How could they work together to produce power? We'll be looking at these questions and more as we learn about wind and hydrogen fuel cell power.

## Method

---

You will be using your miniature turbine and a fan to simulate a full-sized wind turbine in a consistently windy location. You will change the characteristics of your turbine to make it more effective at generating electricity and compare its power output to a hydrogen fuel cell.

## Procedure

---

Use these guiding questions as you assemble and test your turbine. Record anything interesting you observe while building and experimenting in the **Observations** section below.

1. Look at the three different types of blades available (labeled A, B, and C). How are they similar? How are they different? Discuss with your group which type of blade you think would work best with your turbine.
2. Select the type and number of blades you want to test. Why do you want to test this type of blade first? Do you think it will be better or worse than the other types?
3. Check that the blades are in the same position using the three notches near the white bases of the blades. Rotate the individual blades if needed to get all the blades into the same position. Would your turbine still work if the blades were in different positions?
4. Insert the blades into the Rotor Base and put the Blade Holder and the Blade Assembly Lock, then

attach the Blade Unit to the metal shaft of the turbine. Can your blades be positioned backwards?

How do you know if there's a "right way" for a blade to be positioned?

5. Connect the base of the turbine to the LED lights using the black and red wires. Why do you think the lights need two wires to work?
6. Turn on the fan and position it in front of the turbine. It will work best if you keep the fan close to the turbine and line up the center of the fan with the center of the turbine. Why would changing the position of the fan affect the wind hitting the turbine?
7. Record your observations in the data table on the next page: Did the lights turn on? Were they dim or bright? .
8. Discuss what you observed with your group and discuss what you want to change: the number of blades, the angle of the blades, the type of blades, or some combination of those.
9. Repeat steps 1-8 with as many changes as you can think of.

**Observations:**

**Data Table:**

<b>Blade Type (A,B,C)</b>	<b>Number of Blades</b>	<b>Blade Angle (6°, 28°, 50°)</b>	<b>Observations:</b>

**Experimentation: Hydrogen Fuel Cell**

You’ve used the turbine to generate electricity and provide power to lights, but you can also use it to generate hydrogen gas using the hydrogen fuel cell. Disconnect the lights from the base of the turbine and connect the turbine to the fuel cell.

The more electricity the turbine generates, the faster the cylinders will fill with hydrogen and oxygen gas. Use a stopwatch to measure the amount of time it takes to fill the hydrogen cylinder (it will fill up first). Stop the timer when you see bubbles start to escape from inside the cylinder and disconnect the fuel cell. Try different blade characteristics to see which generates electricity fastest. Record your observations on the next page.

Once you’ve generated hydrogen, you can use the fuel cell to create electricity. Attach the motor or the lights to the fuel cell and watch as it puts hydrogen and oxygen back together to make water and create electric current.

<b>Blade Type (A,B,C)</b>	<b>Number of Blades</b>	<b>Blade Angle (6°, 28°, 50°)</b>	<b>Time to fill H<sub>2</sub> cylinder (sec)</b>

Which blade characteristics were most effective at producing electricity? Why do you think it was?

**Experimentation:** Revolutions per Minute

Put a piece of masking tape around one of your blades. This will allow you to count the number of times that blade rotates in a given amount of time. Have one group member operate the stopwatch and another stand next to the turbine to count the revolutions. (It’s easier to count if you stand beside the turbine and count the times the taped blade comes over the top.)

Look at your data from the previous experiment. Note which trials had the brightest lights and which had the dimmest or none at all. Use these extreme cases to find a range for your turbine’s revolutions per minute (RPM).

Have the time keeper say “Go!” once the fan has got the turbine turning and count the number of times the taped blade completes a revolution. Stop counting at exactly 15 seconds. You can then multiply your count by 4 to get the total revolutions per minute. Record your data on the next page.

<b>Blade Type (A,B,C)</b>	<b>Number of Blades</b>	<b>Blade Angle (6°, 28°, 50°)</b>	<b>Count (15 seconds)</b>	<b>Revolutions per minute (count × 4)</b>

What was your greatest number of revolutions per minute? How did that compare with other groups in your class?

**Experimentation:** Fan Distance

Increasing the distance between the fan and the turbine will decrease the energy of the wind that reaches it. What is the maximum distance the fan can be moved and still light up the LEDs?

Look at your data from your previous experiments. What are the trials that lit up the LEDs the brightest? Pick the three or four brightest to use in this experiment.

Start with the fan very close to the turbine and gradually move it backwards until the lights no longer light up. Measure the distance on your desk from the turbine blades to the fan in cm and record your data on the next page.

<b>Blade Type (A,B,C)</b>	<b>Number of Blades</b>	<b>Blade Angle (6°, 28°, 50°)</b>	<b>Maximum Distance (cm)</b>

Which type, number, and angle of blades could keep the lights on for the farthest distance? What do you think that means?

## Analysis

Make a *scientific claim* about the turbine blades: What characteristic was most important in getting them to turn faster? To help you write a claim statement, see [Stating a Scientific Claim](#).

### Claim

What *evidence* can you use from your **observations** and your **experimentation** to back up your claim?

### Evidence

State the *reasoning* you used to make your claim.

### Reasoning

Did other people come to a different conclusion about their turbine blades? Use your observations to *develop an argument* that supports your position and defend your findings.

## Conclusion

---

You were able to change the angle, shape, and number of your turbine blades. *Design an experiment* that would change a feature of your turbine that you weren't able to test.

Based on your data, what is the best angle, shape, and number of blades for your wind turbine?

*Design an experiment* that tests whether different types, angles, and numbers of blades work better in different wind speeds.

How could you use a hydrogen fuel cell with a wind turbine to generate power even when there isn't any wind?

What characteristics do you think would make a good location for wind energy? Do you think your school is located in an area that could provide wind energy? Why or why not?

Based on your observations, do you think there are more benefits of wind power than drawbacks? *Develop an argument* to defend your position and discuss your conclusion with your classmates.

Could you generate more power with more turbines? Explain how you would connect multiple turbines together if you wanted to create more electric current.

Do you think a hydrogen fuel cell is better or worse for generating electricity than a wind turbine? *Develop an argument* to defend your position and discuss your conclusion with your classmates.

## Measurement

---

This section requires a multimeter or the Horizon Renewable Energy Monitor. For help setting up a multimeter, see [Using a Multimeter](#).

If the turbine is spinning faster or slower, we should be able to figure out how much faster or slower it's moving based on how much electricity it's producing. Using a multimeter like the Horizon Renewable Energy Monitor we can measure the voltage and amperage of the electric current. See [Measuring Current in a Circuit](#) if you need to know the difference between volts and amps.

With your multimeter attached to the turbine and LEDs, you should be able to see differences in the current produced by the turbine. Change the number, type, or angle of the turbine blades and see what changes you can observe when the turbine is placed in front of the fan. Record your observations on the table on the next page:

<b>Blade Type (A,B,C)</b>	<b>Number of Blades</b>	<b>Blade Angle (6°, 28°, 50°)</b>	<b>Volts</b>	<b>Amps</b>